



## Product - Information



## Kaps 900

### Cost-effective entry level Dental Microscope with integrated LED illumination

The Kaps 900 is a high performance and user-friendly dental microscope. It supported the dentist during his daily work.

#### Features at a glance

- 5-step magnification changer
- 0-210° inclinable binocular tube  $f = 182$  mm
- Objective  $f = 250$  mm with fine focusing
- Widefield eyepieces  $12.5 \times V$
- integrated LED illumination with fibre optic cable (life time app. 50,000 hours)
- compact mobile stand, h-style
- ergonomic handle

Kaps offers a wide range of **accessories**:

- Kaps monolateral beam splitter with rotation device
- Kaps HD photo adapter for SONY NEX-5/-6/-7 and Alpha 5000/6000
- Kaps Varilex-100 objective  $f = 220 - 320$  mm
- Lens protection
- Sterilization cover for ergonomic handle

