




Product - Information



Kaps TFT-display mount

TFT-display mount for Kaps SOM 62 or mobile stand

With the Kaps TFT-display mount the user can be easily position the display for comfortable viewing or even move it out of the way when not in use.

- ideal ergonomic comfort for the user
- easy and simple mounting on SOM 62 or mobile stand
- adjustable in height above the column
- swing, tilt or turn the flat panel monitor for maximum efficiency
- VESA® Flat Display Mounting Interface (FDMI™) Standard compliant for easy monitor attachment (up to 22“)
- including mounting-set and instruction
- ➔ order no 162.205-000



Displays

Neovo DR-17 P

durable, versatile, high performance display designed expressly for medical environments

- medical certification according to EN60601-1 und IEC60601-1

- 17"-display
- 1280 x 1024 resolution
- NeoV™ optical glass
- AIP technology: PIP; built-in image enhancer
- durable metal casing
- S-Video, CVBS, analogue, digital inputs
- built-in speakers (1.5W x 2)
- DC 24V medical grade power adaptor allows longer distance connection

➔ order no 711-147



Neovo X-17 P

brilliant display that delivers advanced design, durability and versatile performance for progressive professional environments

- 17"-display
- 1280 x 1024 resolution
- NeoV™ optical glass
- AIP technology: PIP; built-in image enhancer
- durable metal casing
- S-Video and CVBS inputs
- built-in speakers (1.5W x 2)

➔ order no 711-164



Neovo X-22

given its hard glass protection and metal casing, the X-22-W promises to excel in all sorts of public environments, however demanding they may be

- 22"-display
- LED-backlit technology with high 1920 x 1080 resolution (full HD)
- NeoV™ optical glass
- Versatile inputs: VGA, DVI, HDMI, S-Video, CVBS (RCA x 2), Audio in/out
- AIP technology: PIP and PBP functions; 3D Comb Filter/Deinterlace/Noise Reduction; image rotation function
- Durable metal housing
- built-in speakers (2W x 2)

➔ order no 800-000-020

